




## LAP TOP CUSHIONS



**A**  
**FULL-ARM WITH CUTOUTS**



**B**  
**FULL-ARM NO CUTOUTS**



**C**  
**HALF-ARM**

### FOR FULL ARM WHEELCHAIRS

Description	Thickness	Chair Width	Reorder #	Black	Mauve
With Cutouts	4"	16" to 18"		<input type="checkbox"/> 307010	<input type="checkbox"/> 307011
With Cutouts	2½"	16" to 18"		<input type="checkbox"/> 307012	<input type="checkbox"/> 307013
With Cutouts	3¼"	16" to 18"		<input type="checkbox"/> 307014	<input type="checkbox"/> 307015
No Cutouts	4"	16" to 18"		<input type="checkbox"/> 307020	



### FOR HALF ARM WHEELCHAIRS

Thickness	Chair Width	Reorder# Black	Mauve	Chair Width	Reorder# Black	Mauve
4	16"	<input type="checkbox"/> 307030	<input type="checkbox"/> 307031	18"	<input type="checkbox"/> 307040	<input type="checkbox"/> 307041
2½"	16"	<input type="checkbox"/> 307032	<input type="checkbox"/> 307033	18"	<input type="checkbox"/> 307042	<input type="checkbox"/> 307043
3¼"	16"	<input type="checkbox"/> 307034	<input type="checkbox"/> 307035	18"	<input type="checkbox"/> 307044	<input type="checkbox"/> 307045

### PURPOSE

Lap Top Cushions are intended to serve as reminders for residents to ask for assistance before leaving the wheelchair. They also help the forward-leaning resident maintain an upright sitting posture by providing a support foundation for his/her arms or elbows. Lap Top Cushions provide a convenient surface for the resident to place personal items or rest his/her hands while holding reading material. **Important notice: Both Full-Arm and Half-Arm Lap Top Cushions with cutouts may be used to control sliding. Do not use the Lap Top without cutouts (#307020) for this purpose.**

### NOTICE

A Lap Top Cushion generally will not be considered a physical restraint if your resident can demonstrate that he/she can remove it at all times. The Lap Cushion is considered to be a restraint for the resident who cannot remove it.

### APPLICATION

1. Lock wheelchair brakes and rotate footrests to the side of the wheelchair.
2. Assist resident into the wheelchair and seat him/her as far back on the seat as possible.
3. Position the footrests to the front of the chair and place the resident's feet on them.
- 4A. (Full-Arm Lap Top With Cutouts) Slide Lap Cushion under the armrest support bar and position the cutouts as shown in illustration A above.
- 4B. (Full-Arm Lap Top Without Cutouts) Slide Lap Top Cushion under the armrest support bars making certain the straight side faces away from the resident and the contoured side is positioned at the resident's waist as shown in illustration B above.

**4C. (Half-Arm Lap Top)** Place Half-Arm Lap Top Cushion between armrest support bars with zipper-side facing away from resident. Secure by positioning cutouts as shown in illustration C above.

5. If resident is unable to remove the Lap Top Cushion and ambulate on his/her own, the resident should be released at least every two hours for exercise, toileting, and repositioning.

**IMPORTANT**

Lap Top Cushions are intended only as aides for resident positioning and safety. They are not substitutes for proper supervision. Residents should be checked frequently and the Lap Top Cushion readjusted if necessary.

**CLEANING**

Lap Top Cushions may be cleaned with a spray-type cleanser and a damp cloth or sponge. A mild disinfectant may be used if necessary. Do not launder cover.