


Bariatric Gel-Foam Cushion



Reorder #	Size (inches)	Reorder #	Size (inches)
<input type="checkbox"/> 915110	20 x 18 x 3 ½	<input type="checkbox"/> 915130	24 x 18 x 3 ½
<input type="checkbox"/> 915120	22 x 18 x 3 ½	<input type="checkbox"/> 915135	26 x 18 x 3 ½

Effective Design

The two-chamber Bariatric Gel-Foam Cushion combines the comfort and therapeutic effectiveness of gel and foam. The high-density foam base prevents bottoming-out, provides stability, and increases comfort. The water-based gel cools the skin to prevent perspiration while relieving interface pressure.

New Low-Shear Cover

The outer cover is made of two-way stretch Low-Shear fabric. It moves with the wheelchair user to reduce skin-damaging shear and friction. The Low-Shear cover is completely launderable.

Important Safety Feature

The Low-Shear cover has a safety strap that prevents the cushion from slipping off the seat. This strap must be secure to the wheelchair frame.

Easy Care

The Low-Shear cover may be washed and dried at temperatures that do not exceed 180°F (82°C). Do not use bleach. The incontinent-proof vinyl inner cover is factory sealed to prevent urine from soaking through to the foam. Spray-and-wipe cleanser and mild disinfectant may be used to clean the vinyl.