

# Contour Cushion™

**Low-shear cover** prevents skin-damaging friction. Incontinent-proof. Removes for easy laundering.

**Safety straps** secure cushion to wheelchair.

**Visco-elastic foam pad** provides support and pressure relief beneath the coccyx.

**Elevated sides** improve lateral stability.

**Ample high density foam** relieves pressure and increases comfort.

**Raised front** controls forward thrust and maintains pelvic neutrality to reduce sacral sitting problems.



**Available gel pod** provides additional pressure relief in the critical coccyx area.

Reorder #	Description	Size	Unit
<input type="checkbox"/> 751810	Contour Cushion™ Foam Only	18" x 16"	Each
<input type="checkbox"/> 751820	Contour Cushion™ Gel-Foam	18" x 16"	Each
<input type="checkbox"/> 751830	Contour Cushion™ Foam Only	16" x 16"	Each
<input type="checkbox"/> 751840	Contour Cushion™ Gel-Foam	16" x 16"	Each

## Purpose

Designed for positioning and pressure relief. The raised front controls forward thrust to maintain pelvic neutrality. Elevated sides promote lateral stability. Ample foam cushioning relieves pressure and adds comfort.

## Easy Care

The cloth cover may be wiped clean using a mild spray cleanser and damp cloth, or it may be removed for laundering. Wash and dry at temperatures below 180°F (82°C). Do not use bleach.