


**Skil-Care™**29 Wells Avenue  
Yonkers, NY 10701914 963-2040  
800 431-2972

## GRIP-EZE

### FINGER CONTRACTURE CUSHION

Re-Order #	Unit
☐ 201140	6/Pk.
☐ 201150	36/Cs.

The GRIP-EZE hand cushion is designed to reduce between-the-fingers irritation, control hand contracture, and prevent patients from digging their fingernails into their palms. The individual finger loops reduce skin damage that results when fingers are rubbed together or forced together by contraction. The resilient sheepskin stuffing allows a mild form of hand exercise.

**LAUNDERING INSTRUCTIONS** Completely launderable. The GRIP-EZE may be washed and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

**Skil-Care Corporation . . . Advanced design for nursing and therapy.**