


HEELS-OFF

Reorder #	Size	Unit
503150	24" x 16" x 7" x 5½"	Each



Purpose

Heels-Off is designed to off-load both heels to prevent pressure sores from forming or aid in the healing of existing ulcers. It elevates the legs to relax spine and prevent leg and lower back pain.

Easy to Use

Place the higher end (7 inches) of the Heels-Off device under the patient's knees. The heels should be suspended from the lower (5½ inches) end.

Easy Care

The low-shear cover wipes clean with a mild cleaning solution and non-chlorine disinfectant. Do not launder.

