



LEG SEPARATOR PADS

Re-Order #	Size	Unit
<input type="checkbox"/> 503420	9" x 15" (Thigh)	Each



Re-Order #	Size	Unit
<input type="checkbox"/> 503460	15" x 15" (Knee)	Each



Skil-Care's Leg Separator pads eliminate the skin irritation and soreness that often appear between your patient's thighs or knees. By placing a soft barrier of plush synthetic sheepskin between knees and thighs, these pads reduce tissue-damaging friction and moisture build-up. Effective in the bed and wheelchair, the pads stay securely in place with non-slip buckles. For complete leg protection, use the Full Leg Pad (Re-Order #503480).

LAUNDERING The Leg Separator Pads may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

Skil-Care Corporation . . . When Safety Matters and Comfort Counts.