

Visit us online
at Skil-Care.com

Skil-Care[™]
CORPORATION

29 Wells Avenue
Yonkers, NY 10701

914-963-2040
800-431-2972



Skil-Care[™]
CORPORATION

POUCH ARM SLING

<u>Re-Order #</u>	<u>Size</u>
<input type="checkbox"/> 902011	Small-Medium
<input type="checkbox"/> 902012	Medium-Large

Skil-Care's durable, pouch-style sling provides uniform weight distribution for either arm. The completely adjustable strap features a secure, non-slip buckle.

LAUNDERING INSTRUCTIONS

Skil-Care's Pouch Arm Sling may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.

Advanced design for nursing and therapy.[™]