

Visit us online  
at Skil-Care.com

**Skil-Care**<sup>™</sup>  
CORPORATION

29 Wells Avenue  
Yonkers, NY 10701

914-963-2040  
800-431-2972

# SLIP-GRIP

Re-Order #	Color	Size	Unit
<input type="checkbox"/> 909250	Blue	12" x 120"	Roll
<input type="checkbox"/> 909251	Green	12" x 120"	Roll
<input type="checkbox"/> 909252	Raspberry	12" x 120"	Roll

## Purpose

Slip-Grip is a nonslip matting that may be used to prevent cushions from slipping off chairs; hold plates so that unsteady eaters can feed themselves; and position instruments on uneven surfaces. Sold in 12: x 120" rolls, Slip-Grip may be cut to required size.

## Important Safety Information

Slip-Grip loses its nonslip quality on wet surfaces. Do not use Slip-Grip on shower or bathroom floors. To avoid skin shearing, do not place Slip-Grip in direct contact with resident's skin.

## Cleaning

Slip-Grip may be washed with soap and water and air dried. Do not machine wash and dry. Do not autoclave. Make certain that Slip-Grip is completely dry before using. Slip-Grip can EO sterilized.

**Skil-Care**<sup>™</sup>  
CORPORATION

*Advanced design for nursing and therapy.*<sup>™</sup>