


SOFT BOOT

Re-Order #	Previous #	Size	Unit
□503160	8005	Universal	Pair



Skil-Care's Soft Boots offer extended decubitus protection that provides pressure relief for your patient's heel, ankle, and lower calf. Made from plush synthetic sheepskin, these boots are comfortably secured with hook-and-loop closures. Large vent holes promote air circulation to reduce perspiration and moisture build-up.

LAUNDERING INSTRUCTIONS

The Soft Boots may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended. To prolong the useful life of the hook-and-loop closures, secure them before laundering. This will prevent the hook material from picking up lint during the washing and drying process.