


CONFORM Thin-Line™ CUSHION

| Reorder # | Description | Size |
|-----------|-----------------|-----------------|
| 753150 | ConForm Cushion | 18" x 16" x 2¼" |

Skil-Care's Thin-Line™ cushions are low-rise cushions that enable many shorter users to reach the floor or wheelchair footrests.

Purpose

The special temperature-sensitive foam used in the ConForm Cushion molds to the dimensions to provide uniform pressure relief. It is designed to promote better sitting posture by reducing forward sliding.

- 
- **Safety Strap** secures cushion to wheelchair
 - **Low-Shear II cover** prevents skin-damaging friction. Fluid-proof. Wipes clean.
 - **Low-rise barrier** prevents sliding by controlling forward thrust and maintaining pelvic neutrality.
 - **Contour Leg troughs** promote thigh alignment for postural symmetry.
 - **Visco-elastic foam** provides uniform pressure relief by conforming to user's dimensions. High and low spots are eliminated.

Easy Care

The Low-Shear II cover wipes clean using a damp cloth and mild cleanser. Do not launder.