


Ultra-Soft Heel Cushion

Reorder	Size	Unit
503030	Universal	Pair



Protects sensitive heels against pressure sores. Filled with fluffy spun fiber for cushioning and comfort. Covered with low-friction brushed polyester fabric to reduce skin shear. An extra-wide strap with Velcro® closures keeps the Ultra-Soft Heel Cushion securely in place.

Application

For maximum pressure relief, do not allow heel to rest on seam. Place Ultra-Soft Heel Cushion on foot according to instructions on product label.

Laundering

The Ultra-Soft Heel Cushion may be washed and dried at temperatures that do not exceed 180°F (82°C)