



## WHEELCHAIR PLUSH PADS (For Seat and Backrest)

Re-Order #	Previous #	Unit
703050	9070	2 Per Set

These pads add a touch of luxury to any wheelchair. The high-pile synthetic sheepskin provides pressure relief plus improved air circulation. Packaged as a set of two units, one for the seat and one for the backrest, each pad has tie straps that secure it to the wheelchair.

### LAUNDERING INSTRUCTIONS

SkiL-Care's Plush Pads may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended.