

# Lateral leaning causes ischial pressure and uncomfortable positioning. Skil-Care's Adjustable Lateral Support corrects side-leaning for proper pressure redistribution and greater comfort.



*“Pressure ulcers are usually located over a bony prominence, such as the ischial tuberosity.” (F14) Side-leaning wheelchair users place increased pressure on the vulnerable ischials. The Adjustable Lateral Support keeps residents upright to maximize weight distribution and minimize pressure ulcer risk.*



## The Adjustable Lateral Support includes these important features:



Complete adjustability to accommodate a range of body widths.



Comfortable foam-cushioned side pads that control lateral leaning. The pads are height adjustable to stay below the axilla and not interfere with user's ability to reach hand-rims on wheels.



Cushioned backrest that adds comfort and back support.

Reorder #	Size	Width Adjustment	Unit
706050	Small	8" - 12"	Each
706055	Medium	12" - 18"	Each
706059	Large/Bariatric	18" - 24"	Each

**Skil-Care's Adjustable Lateral Support is available in three sizes for petite to bariatric residents. It is adaptable to all wheelchairs and wipes clean for easy care.**

**For Information Contact:**

**Skil-Care™**  
CORPORATION

29 Wells Ave., Yonkers, NY 10701

Advanced design for nursing and therapy.™

914-963-2040 800-431-2972 Fax 914-963-2567