

Heel-Float™ Plus+ for 5-point pressure reduction.

Advanced design makes Skil-Care's new Heel-Float™ Plus+ the only heel suspension boot that reduces or eliminates pressure in these five critical areas...



1. **Toes** Extended foot support keeps sheets and blankets off toes.
2. **Bunion** Cloth-covered Velcro® closure secures below bunion to avoid pressure and the pain it causes.
3. **Malleoli** Cutouts for the lateral and medial malleoli provide pressure relief and eliminate skin shear.
4. **Heel** Suspension design off-loads heel over air cavity for zero pressure.
5. **Ankle** Water-based gel pack reduces pressure under ankle and cools to prevent skin-damaging perspiration.

Heel-Float™ Plus+ provides these important features...

Air-Flow Design



Prevents perspiration and aids in healing existing sores.

Rigid Insert



Provides support to prevent foot drop.

Cozy Cloth™ Cover



Low-friction cover prevents skin shear and wicks away moisture. **Cozy Cloth™** is sewn to Velcro® closures to prevent skin abrasion.

Skil-Care's Heel-Float Plus+ is completely launderable for maximum infection control and is replacement guaranteed for six months.

Reorder #	Size	Unit
503044	Small (3" Wide)	Each
503045	Medium (4" Wide)	Each
503046	Large/Bariatric (5" Wide)	Each

Skil-Care™
CORPORATION

29 Wells Ave., Yonkers, NY 10701

Advanced design for nursing and therapy.™

914-963-2040 800-431-2972 Fax 914-963-2567

For Information Contact: