

Heel-Float™ off-loads the heel for total pressure relief.

Skil-Care's Heel-Float™ suspends the heel over an air cavity to prevent pressure sores or aid in the healing of existing ulcers. Make this zero-pressure suspension boot part of your F314 care plans.

**Reusable
Launderable
Guaranteed for
6 months**



- New Cozy Cloth™** cover comfortably reduces skin shear and wicks away moisture.
- Anti-foot-drop** construction controls dorsal flexion.
- Water-based gel insert** reduces pressure beneath the ankle and cools to reduce perspiration. (May also be used for hot/cold therapy.)
- Foam cushioned sides** protect the lateral and medial malleoli for side-lying patients.
- Air cavity** floats the heel for zero pressure. Provides air circulation to aid healing of existing ulcers.

Available in three sizes for better fit and effectiveness.



No ankle or shin abrasion
Heel-Float™ protects sensitive skin by placing a **Cozy Cloth**™ barrier between the Velcro® closures and the patient's ankle and shin.



- Foot Stabilizer**
- Use with Heel-Float to prevent foot rotation
 - Recommended for hip-fracture patients
 - Durable and reusable
 - Fits all makes of heel-suspension boots

Reorder #	Description	Size	Unit
503034	Heel-Float	Small (children and petite adults)	Each
503035	Heel-Float	Medium (average adult)	Each
503036	Heel-Float	Large (large adults and bariatric patients)	Each
503039	Foot Stabilizer	Small/Medium	Each
503038	Foot Stabilizer	Large/Bariatric	Each

All models of the Heel-Float and Foot Stabilizer fit both the left and right foot.

For Information Contact:

Skil-Care™
CORPORATION

29 Wells Ave., Yonkers, NY 10701

Advanced design for nursing and therapy.™

914-963-2040 800-431-2972 Fax 914-963-2567