

# Hip-Ease™ – More cushioning where it counts.



**Velcro™ waist and leg bands** for easy-on/easy-off convenience.

**Air-Flo fabric** that's cool, comfortable and completely launderable.

**High-density foam pads** are securely positioned over the vulnerable upper femur.

**Ample, impact-absorbing foam pads decrease fall-related hip trauma.**

**Skil-Care's new Hip-Ease™ ...**

- ▶ Fits over all undergarments, including incontinence briefs
- ▶ Features removable/replaceable, incontinent-proof, vinyl-covered foam pads
- ▶ Comes in six sizes for better fit and full coverage of hip area

Reorder #	Size	Waist Measurement		Unit
911450	X-Small	26" - 28"	66 - 71 cm	Each
911452	Small	28" - 30"	71 - 76 cm	Each
911454	Medium	30" - 34"	76 - 86 cm	Each
911456	Large	34" - 38"	86 - 96 cm	Each
911458	X-Large	38" - 42"	96 - 106 cm	Each
911460	XX-Large	42" - 46"	106 - 116 cm	Each
911470	Replacement Foam Pads Vinyl Covered	_____	_____	Pair

For Information Contact:

**Skil-Care™**  
CORPORATION

29 Wells Ave., Yonkers, NY 10701

Advanced design for nursing and therapy.™

914-963-2040 800-431-2972 Fax 914-963-2567