

• Impact Absorption • Step-on Stability

SafeSide Fall Mats provide both to protect your residents against fall-related injuries.



Skil-Care's **SafeSide Fall Mat** uses shock-absorbing bonded foam to cushion falls from bed. Bonded foam resists compression under impact to provide greater protection against fall-related trauma.

The **SafeSide Fall Mat** is designed with step-on stability to prevent tumbles when a resident walks on or approaches the mat. A slip-resistant vinyl bottom keeps the SafeSide Mat securely in place. The dense bonded foam won't "bottom-out" when stepped on –this means safer footing for residents. The low-profile/tapered-edge design reduces the likelihood of a resident or caregiver stumbling as he or she approaches the mat. And the tapered edge facilitates wheelchair roll-on.



The durable vinyl covers of **SafeSide Fall Mats** have heat-sealed seams to protect the foam against moisture and bacteria –an important infection control feature.

SafeSide Fall Mats are available in your choice of 24-inch- and 36-inch-wide bi-fold or non-folding models. All mats are 68 inches long and fit easily under all hospital-style beds.

Reorder #	Style	Size	Unit
911570	Non-Folding	68" x 24" x 1"	Each
911572	Bi-Fold	68" x 24" x 1"	Each
911575	Non-Folding	68" x 36" x 1"	Each
911577	Bi-Fold	68" x 36" x 1"	Each

Available option. A special strap that permits the bi-fold models to be secured to the bed footboard or hung from a wall hook is available. The order number for the strap is 911539.

For Information Contact:

Skil-Care[™]
CORPORATION

29 Wells Ave., Yonkers, NY 10701

Advanced design for nursing and therapy.[™]

914-963-2040 800-431-2972 Fax 914-963-2567